

# SLEEK AND SENSIBLE. SIMPLIFIED.

The new, sleek X-Act ROM™ Lite is your sensible solution to post-op knee bracing. This lightweight, easy-to-use brace features the durable X-Act ROM hinge, full circumferential foam, and an easy-to-adjust strap system. The simplified telescoping uprights can fit a variety of patients and can be contoured to the patients' leg to provide varus and valgus support.





#### EASY TO USE STRAPPING SYSTEM

Velcro strapping for quick and easy application. Strap extension kit provided to fit a variety of patients.

### ADJUSTABLE UPRIGHTS

Easy Push Button System has telescoping adjustability between 18in-25.5in (46cm-65cm) can also be bent to provide varus/valgus support.

#### **FULL CIRCUMFERENTIAL FOAM**

Cool, breathable foam provides a comfortable fit to patients' leg and helps minimize migration of brace.



The X-Act ROM is DonJoy's lightest post-op brace, allowing a comfortable, no-slip fit with precise protection during healing.



## **Ordering Info**


Part # 11-2161-9	X-Act ROM Knee Lite, Universal
Part # 25-6755	X-Act ROM Knee Hinge Zip Ties

Call 800.336.6569 or visit www.donjoy.com



DJO, LLC | A DJO Global Company

**T** 800.336.6569 **D** 760.727.1280 **F** 800.936.6569 1430 Decision Street I Vista, CA 92081-8553 I U.S.A. DJOglobal.com



DJO Global proudly sponsors:









